

Jazzy[®]

SELECT 14 SERIES

THE EASIEST POWER CHAIR TO USE.
#1 IN MANEUVERABILITY.

The Jazzy[®] Select[™] 14 Series delivers the advantages of patented mid-wheel drive for maneuvering in tight, compact spaces. Both the Jazzy Select 14 and Jazzy Select 14 XL are loaded with performance and convenience features like a powerful 60 amp PG VR2 controller, larger batteries for greater per-charge range, and on-board battery charger.



Jazzy[®] Select[™] 14

- 300 lbs. weight capacity
- Maximum contoured dimensions of up to 20" x 20"

#1 Physician Approved



Jazzy[®] Select[™] 14 XL

- 400 lbs. weight capacity
- Contoured seating available up to 22"W x 22"D



The Jazzy® Select™ 14 series models include:

- Jazzy Select 14-C
- Jazzy Select 14-SS
- Jazzy Select 14 XL-C
- Jazzy Select 14 XL-SS

FEATURES

- Patented mid-wheel drive design for tight space maneuverability
- 14" knobby tires and optimum weight distribution provide superior traction
- NF-22 batteries for maximum range
- New 60 amp, PG VR2 controller
- On-board, 4 amp charger

OPTIONS

- Contoured high-back seat ¹
- Contoured high-back with solid seat pan ¹
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Safety flag

¹ Standard on Jazzy Select 14 XL

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %.

	Jazzy® Select™ 14	Jazzy® Select 14™ XL
Weight capacity:	300 lbs.	400 lbs.
Maximum speed:	Up to 4.5 mph	Up to 4.5 mph
Ground clearance:	3"	3"
Turning radius:	21"	21"
Overall length:	38" without front riggings	38" without front riggings
Overall width:	24.75"	24.75"
Wheels:	Drive: 14"	Drive: 14"
Front anti-tips:	6" solid	6" solid
Rear casters:	8" solid	8" solid
Suspension:	Limited	
Drive train:	Two-motor, mid-wheel drive	
Braking system:	Intelligent braking (electronic, regenerative disc brakes)	
Std. Electronics:	PG VR2, 60 amp	
Specialty controls:	N/A	
Battery charger:	On-board, 4 amp	
Per-charge range:	Up to 15 miles	
Battery requirements:	(2) 12 volt, deep cycle	
Battery size:	NF-22	
Transport weight		
Battery weight:	38 lbs. each	38 lbs. each
Base weight:	105 lbs.	106.5 lbs.
Std. Seat dimensions:	20" x 20"	22" x 22"
Warranty:	5-year limited warranty on frame; 1-year limited warranty on electronics; 1-year limited warranty on drive motors	

Colors:



Pride SilverStar Lifts & Ramps

- Full Platform
- Backpacker® Plus
- Power chair lift with straps
- Silver Boom 300

Full Platform



Backpacker® Plus



These lifts are compatible with the Jazzy® Select™ 14 Series.

Pride[®]
Mobility Products Corp.