

Backpacker™ PLUS

All-black battery box, bracket, and new matte finish



One-piece, fully assembled drive-on/drive-off platform



Single-piece horizontal drive rail: Minimum welding points for noise reduction, increased rigidity, and more fluid movement



Removable, front drive panel for easier serviceability



New, angled, low-profile platform



4" more usable space

New, adjustable wheel chocks on the platform surface



Easy to install and use. The new Pride® Silver Star™ Backpacker™ Plus Interior Platform Lift features a pre-assembled single-piece platform that you can drive your scooter or power chair on or off from either direction. Enhanced construction and a fresh new look give the Backpacker Plus the quality you expect and deserve from Pride Silver Star Lifts.

- Designed and built by Pennsylvania workers
- The quality construction and design you expect from Pennsylvania-made Pride® Mobility Products
- 325 lbs. lifting capacity
- Easy-to-use 2-button hand control
- 3-year transferable warranty
- Can lift a variety of Pride® Scooters and Power Chairs

Pride® SilverStar™ Lifts

182 Susquehanna Ave. • Exeter, PA 18643 • www.silverstarmobility.com • www.pridemobility.com
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)