

Full Platform Exterior Lift

Shock-absorbing "cradle effect" protects power chairs from jarring vibrations.



New drive-on/ drive-off platform



Easy to use, three-point power chair tie down system allows for one handed operation and is self-retracting and self-tensioning.



All new design dramatically reduces metal to metal contact for enhanced performance and longer product life.

Innovative, lightweight self-leveling platform allows for easy loading and unloading, even on inclines or uneven terrain.

Your on-the-go lifestyle demands that you bring your power chair with you. Now, Pride Mobility Products® Corp., the number one power chair and scooter manufacturer, offers the lift you need to take your mobility product on the road. Backed by the service you have come to expect and trust from Pride®, with the full support of our expert technical service staff, a Pride® Silver Star™ Lift is the natural complement to your Pride Power Chair.

A Pride Silver Star exterior lift makes it easy to take your Pride Power Chair with you wherever you go, complementing your product perfectly and enhancing your mobility. An all steel frame and protective coating mean Pride Silver Star Lifts are built to be durable, while quality construction means you can count on your lift to be dependable. From loading to lifting, Pride Silver Star Lifts are designed with ease of use in mind.

- 325 lbs. lifting capacity
- 3-year transferable warranty
- Most adaptable exterior lift available

NEW! • Drive-on/drive-off platform allows user to enter from either direction

Pride® SilverStar™ Lifts

182 Susquehanna Ave. • Exeter, PA 18643 • www.silverstarmobility.com • www.pridemobility.com
(800) 800-8586 (U.S.) • (888) 570-1113 (Canada)

Full Platform Exterior Lift



Unoccupied lift stowed.



Lift down, ready for product.



Loaded product secured and ready for lifting.



Product lifted and ready for transport.

Base weight lift	92 lbs. ¹
Power chair overall weight	325 lbs.
Overall platform length	49.25" (41" usable surface)
Maximum product width	27"
Maximum product wheelbase	35"
Maximum travel distance from the ground	20"
Motor	Sealed 12 volt DC
Hitch class	Class II or III ^{2,3}
Vehicle platform.....	Full-size van, Mini van, Truck, Sedan, SUV and Crossover
Warranty.....	3-year transferable



Can easily lift the following Pride® Power Chairs ^{3,4}

- Jazzy® 600
- Jazzy® 610
- Jazzy® 1103 Ultra
- Jazzy® 1113 ATS
- Jazzy® Select™
- Jazzy® Select™ 6
- Jazzy® Select™ 7
- Jazzy® Select™ 14
- Jazzy® Select™ 14 XL
- Jet™ 3 Ultra

Can easily lift the following Quantum® Power Chairs ^{3,4}

- Quantum® 600
- Quantum® 610
- Quantum® 1103 Ultra
- Quantum® 1107
- Quantum® 1121
- Quantum® 1122
- Quantum® 6000
- Quantum® 6000 XL
- Quantum® 6000 Z
- Quantum® Dynamo 1107

Options

- License plate kit
- Optional battery pack with charger
- Swing-away joint
- Extended swing-away joint
- Power chair weather cover

- ¹ Weight does not include any accessories
² 190 lbs. maximum product weight with a Class II hitch
³ Vehicle dependent
⁴ Call for details when lifting units with power positioning systems

Notes: The Full Platform Exterior Lift can lift all of the listed power chairs with contoured, Synergy® or Specialty seating. Please call Customer Service when requesting information about lifting a scooter.

Pride® SilverStar™ Lifts 

182 Susquehanna Ave. • Exeter, PA 18643 • www.silverstarmobility.com • www.pridemobility.com
 (800) 800-8586 (U.S.) • (888) 570-1113 (Canada)